

CLEANING INFORMATION

MTC Systems recommends that our products are cleaned regularly according to CDC guidelines and protocols. MTC Systems products are designed for easy cleaning, and to be chemical and spill resistant. Our components are manufactured from high grade, non-porous, chemical resistant materials such as ABS and HDPE plastic, and from smooth, chemical resistant powder-coated steel and aluminum. Work surface components incorporate anti-microbial additives. Although is not recommended to rely solely on those protections for infection control and proper hygiene.

NOTICE: Products are made to withstand repetitive cleaning, but it is important for your facility to spot test the cleaning products to prevent damage. To ensure compatibility, test the cleaning product by applying to an inconspicuous area first, minimizing the time of exposure and the amount of cleaning agent used, following instructions recommended by the supplier, to prevent any damage to the surface. Check for damage to the inconspicuous area before using on other surfaces.

CLEANING METHODS, WIPES, AND SOLUTIONS

Most common disinfectant wipes and diluted cleaning solutions that are safe to use in patient care environments are also acceptable to clean with. Follow manufacturers' guidelines for the specific wipes or solutions being used. Use a soft cloth or wipe to clean the product. To avoid excess dampness and/or damage, do not immerse products in liquid and do not allow liquids to spill inside. Do not use steel wool or other abrasive cleaning methods.

ACCEPTABLE CLEANING

AGENTS INCLUDE:

- Soap and Water
- Diluted Isopropyl Alcohol Solutions (up to 70% alcohol as recommended by the US CDC).
- Diluted Bleach Solutions (carts up to 10%, LCD screens up to 5%)
- Hydrogen Peroxide Solutions
- Quaternary Ammonium Compound (e.g., Sani Wipes or CaviWipes)
- Before disinfection you should remove visible stains and buildup with commonly used cleaners.

UNACCEPTABLE CLEANING

METHODS AND CHEMICALS:

- Non-Diluted Bleach or Chemicals
- Acetone
- Mineral Spirits
- Abrasive Cleaners and Polishes
- Baking Soda
- Paint Thinners or Solvents
- Any Flammable, Toxic, or Corrosive Chemicals
- **Do not use** Industrial Washing Machines above temperatures of 167 °F, which is the level where warping of resins occurs.
- **Do not use** Medical Cleaning Booths